

AB-2970 CLETA-S Quick Manual

Composition (Storage)

CL reagent	1mg-----1 bottle (-20℃)
Substrate solution	20mL-----1 bottle (4℃)
Enzyme solution	8mL-----1 bottle (4℃)
Dilution buffer	8mL-----2 bottles (4℃)

Component

	Main Component
CL reagent	MPEC(powder **1)
Substrate solution	Xanthine
Enzyme solution	Xanthine Oxidase(1.25unit /mL)
Dilution buffer	HEPES

**1 : Please make a stock solution by dissolving the 1mg of MPEC powder with 1.78ml ethanol and divide it into small bottle which for each experience, and store it at -80degreeC.

For working solution, dilute 50μL of stock solution with 150 μL of distilled water.

Summary of Antioxidant Activity Measurment

superoxide is generated from xanthine-xanthine oxidase reaction.

The emission of light is provided by adding a chemiluminescence reagent in generated superoxide.

Function of antioxidant activity is evaluated adding an antioxidant enzyme or an antioxidant material in this reaction.

- * CL reagent + Dilution buffer + Substrate solution

. Negative Control
- * CL reagent + Enzyme solution + Substrate solution

. Positive Control
- * CL reagent + Dilution buffer + Substrate solution + Target sample

. Sample (N)
- * CL reagent + Enzyme solution + Substrate solution + Target sample

. Sample (P)

the antioxidant activity is calculated in below formula.

$\text{Antioxidant activity} = 1 - [(\text{sample(P)} - \text{sample(N)}) / (\text{Positive control} - \text{Negative Control})]$

Operating Method

Recommended volumes for 1 measurement

(1) Sample (or Blank solution)	10 μ L
(2) CL reagent	10 μ L
(3) Enzyme solution (or Dilution buffer)	80 μ L
(4) Substrate solution	200 μ L

{Measurement}

- ① Prepare mixing solution with CL reagent & Dilution buffer.
Mix 10 μ L of CL reagent and 80 μ L of dilution buffer per 1 time of measurement.
The solution should be prepared with over turn mixing in appropriate volume for 1 measurement before using.
(ex: prepare mixing solution with measured sample amount \times four times and additional two times at a time.)
- ② Prepare mixing solution with CL reagent & Enzyme solution.
Mix 10 μ L of CL reagent and 80 μ L of enzyme solution per 1 time of measurement.
The solution should be prepared with over turn mixing in appropriate volume for 1 measurement before using.
(ex: prepare mixing solution with measured sample amount \times four times and additional two times at a time.)
- ③ Inject 10 μ L of blank solution into a measuring container. Prepared and diluted sample is contained in the blank solution.
- ④ Inject 90 μ L of prepared ① mixing solution into ③ the measuring container.
Set the measuring container to the device, and inject 200 μ L of substrate solution.
After measuring of luminescence for 10 sec., the value is considered as [Negative control].
- ⑤ Inject 10 μ L of ③ the blank solution into a new measuring container, and inject 90 μ L of ② mixing solution.
Measure it same as process ④, the value is considered as [Positive control].
- ⑥ Inject 10 μ L of a sample into a new measuring container, and inject 90 μ L of ① mixing solution.
Measure it same as process ④, the value is considered as [Sample(N)].
- ⑦ Inject 10 μ L of ⑥ the sample into a new measuring container, and inject 90 μ L of ② mixing solution.
Measure it same as process ④, the value is considered as [Sample(N)].
- ⑧ Calculate antioxidant activity with the above mentioned formula.

Reference

1. Shimomura, O., Wu, C., Murai, A., and Nakamura, H. (1998)
Evaluation of Imidazopyrazinone-Type Chemiluminescent
Superoxide Probes and Their Application to the Measure-
ment of Superoxide Anion Generated by *Listeria*
Monocytogenes: Anal. Biochem 258, 230-235