


# Transfer Buffer AE-1465 EzFastBlot




- ◎ Simple preparation; Dilute it 10 times with distilled water
- ◎ Methanol-Free  
→ Easy waste liquid treatment
- ◎ Fast electrophoresis transfer
- ◎ Available for Wide range of molecular weight

## Fast electrophoresis transfer



## Available for wide range of molecular weight (~200kDa)



Product	AE-1465 EzFastBlot
Code	2332590
Package	1 Bottle
Volume	500 mL (10x stock solution)
Use	For blotting x 50 times
Main component	Tris, Amino acids
Storage	at room temperature for 1 year



## ATTO Corporation

Providing research solutions for Biochemistry/Molecular Biology/  
 Genetic engineering development  
 Head Office: 3-2-2 Motoasakusa, Taito-ku, Tokyo,  
 111-0041 JAPAN  
 TEL: 81-3-5827-4863 FAX: 81-3-5827-6647  
 Email: eig@atto.co.jp